

WAIVER AND RELEASE OF LIABILITY FOR UV CLEANING & RE-USE OF N95 MASKS

I, the undersigned, am a legal representative of [REDACTED] (“**Entity**”) and hereby state that I have read and understand this Waiver and Release of Liability for UV Cleaning & Re-Use of N95 Masks (“**Release**”) and hereby expressly accept and agree to the following:

Entity understands and accepts that [REDACTED] will provide UV cleaning per the protocol outlined below and as more fully explained in the attached workflow entitled “Preparing Your N95 Masks for UV Cleaning and Re-Use” (“**Workflow**”) hereto and incorporated herein by reference (“**Protocol**”):

- Clean N95 masks will be shipped to the site in white paper bags placed inside sealed white containers and dropped off at the determined locations at each identified site.
- Dirty N95 masks to be placed in provided brown paper bag and placed in sealed in dirty gray container.
- There will be one stop at each site per day to pick up dirty masks.
- All N95 masks should be specific to the individual that uses it and labeled with black sharpie marker.
- Dirty N95 masks will be run through the UV Light Sanitization machine up to five (5) times before the mask is discarded. This will be recorded by tally marks on the masks.
- If an employee at the sterilization facility deems the mask to be unsuitable for use, they can discard it at their discretion.

Entity acknowledges and agrees that it is in the best position to assess its needs, processes and applications with respect to COVID-19 applications. Entity acknowledges that the N95 masks cleaned by [REDACTED] have not been designed for use in healthcare to prevent the spread of diseases or to prevent healthcare personnel from airborne exposure, including COVID-19, which may cause serious or life-threatening disease, including severe respiratory illness. Before using any N95 mask cleaned by [REDACTED] Entity shall bear sole responsibility in assessing, evaluating, determining and validating the N95 masks and whether they are suitable for Entity’s intended use. [REDACTED] makes no warranties and expressly disclaims any warranty relating to the efficacy of any of the N95 masks cleaned by [REDACTED] in preventing the spread and/or contraction of COVID-19.

Entity accepts the N95 masks from [REDACTED] following UV cleaning “**AS IS**”.

Entity further waives, releases, and forever discharges and holds [REDACTED] its successors, predecessors, assigns, officers, directors, employees, agents, representatives and affiliates harmless from any damages, injuries, deaths, losses or claims including attorney's fees and costs incurred related to the possession and use of the UV cleaned N95 masks by Entity, its agents, employees and any EMS workers to whom Entity gives [REDACTED] UV cleaned N95 masks. **ENTITY AGREES AND UNDERSTANDS THAT ITS OBLIGATIONS TO RELEASE, INDEMNIFY AND HOLD [REDACTED] HARMLESS INCLUDES [REDACTED] NEGLIGENCE RELATED TO THE UV CLEANING AND RE-USE OF THE N95 MASKS IN CAUSING IN WHOLE OR IN PART ANY DAMAGES, INJURIES, DEATHS, LOSSES OR CLAIMS.**

Entity understands that this Release limits Entity's and its agents, employees and EMS workers' legal rights associated with the use of UV cleaned N95 masks by [REDACTED] and Entity further understands that it is binding on Entity, its agents, employees, successors and legal representatives. Entity acknowledges that if portions of this Release are deemed invalid by law, Entity agrees that the remaining portions will remain intact and enforceable.

By signing this Release, the undersigned, on behalf of Entity, affirms that the undersigned has read this Release in its entirety, all of the undersigned's questions have been satisfactorily answered, and the undersigned fully understands what he/she has read.

Dated this ____ day of _____, 2020.

[REDACTED]


By: _____
Name:
Title:

Preparing Your N95 Mask for UV Cleaning and Re-use

Please follow these steps below for your safety and to ensure the proper return of your mask.

	Steps	Key Details
1	Label your mask	All N95 masks MUST be labeled to ensure the proper return of your mask, with a black permanent marker.
1a	HS/HS labeling	Label your mask with your site, department, first initial, last name, and date of first use. The location you list will be the location your mask will be returned to once decontaminated.
1b	Other organization labeling	Label your mask with your organization, first initial, last name, and date of first use.
2	Daily use	Please limit the daily donning of new masks as much as possible (extended use is encouraged, see CDC guidelines, https://www.cdc.gov/niosh/topics/hcwcontrols/recommendedguidanceextuse.html). Each health care professional is responsible for ensuring the proper fit and integrity of each mask upon re-use.
3	After use - brown paper bag	All used N95 masks are to be discarded in a brown paper "dirty" bag. Place the bagged mask in your department/organizations "dirty drop off" location in the grey "dirty bin".
3a	Pick up	Couriers will collect the bins from the main pick up location to be brought to the decontamination site on a recurring schedule that will be set based on demand/need.
4	Clean masks returned in white paper bag	Masks sent for decontamination will be returned to you in a new white "clean" bag stapled at the top. All decontaminated N95 masks should be kept in this bag until used.
4a	New brown bag	A new brown bag will be included to be used as your "dirty" discard bag.
4b	Tally marks	Tally marks will be added to the mask by decontamination staff each time the mask undergoes the decontamination process.
5	Reusing a mask	As a general rule, the N95 Masks can be UV cleaned/decontaminated 3-5 times. If an employee at the UV cleaning/decontamination facility deems a mask to be unsuitable for use, they will discard it at their discretion. When donning, make sure to fit the Mask properly and test the breathing capacity.

High-level Process Map



Contact for questions:

Marc Sachtjen, HSHS Manager-Central Services, Central Sterile

Cell: 608-512-9900

Email: Marc.Sachtjen@hshs.org